



AtomoRapid™ HIV

Professional use rapid
diagnostic test

Medical device innovator
Atomo Diagnostics specialises in the design and manufacture of easy-to-use rapid diagnostic test products that have won the respect of frontline healthcare providers worldwide.

AtomoRapid™ HIV is manufactured with strict adherence to the highest international quality standards, including ISO 9001 and ISO 13485, to ensure excellent product performance and reliability.



Excellent Performance

The all-in-one AtomoRapid™ HIV device greatly simplifies the test procedure, reducing user errors and improving test accuracy.




100%
Sensitivity



99.8%
Specificity


All-in-One Convenience

The innovative, user-centred design of AtomoRapid™ HIV removes the need for multiple test components by integrating them into the device platform.



Accurate Blood Collection

The easy-to-use blood collection unit is designed to collect exactly the right sample volume, ensuring optimal test performance.











Built-in Safety Lancet

The auto-retracting safety lancet eliminates the risk of hazardous sharps injuries by locking the needle inside the device after use for safe disposal.

Easy Blood Delivery

Blood is delivered cleanly to the test strip by simply rotating the blood collection unit. The test result is displayed clearly in just 15 minutes.

Quick Reference

-  Low blood volume
-  Result in 15 minutes
-  Suitable for use in POC settings
-  100% Sensitivity
-  99.8% Specificity
-  Store at 2-30°C
-  24 month shelf life
-  Single use only

Small Blood Sample

AtomoRapid™ HIV needs only a small drop of blood to give an accurate result. Smaller blood samples are easier to collect and are more comfortable for the patient.

Results in Minutes

Ideally suited for point of care testing, AtomoRapid™ HIV gives an accurate result in 15 minutes and makes testing simpler, safer and more convenient.

Award Winning Design

The innovative design of AtomoRapid™ HIV has been recognised internationally for its performance and its ease-of-use.



Product code:
MRLF001

atomo diagnostics

SIMPLY BETTER DIAGNOSTICS

Australia

Atomo Diagnostics Pty. Ltd.
Level 2, 701-703 Parramatta Road
Leichhardt, NSW 2040
Australia

T: +61 (0)2 9099 4750

E: info@atomodiagnostics.com

United Kingdom

T: +44 (0)1295 722867

South Africa

T: +27 (0) 31 201 1484

atomodiagnostics.com